
CERTIFICATE OF ANALYSIS

Product	10CFU™ Sensitivity Standard <i>Mycoplasma orale</i>
Cat. No.	102-2003
Lot No.	102o1077
Date of Expiry	31.07.2019
Date of Analysis	18.07.2017
Quantity	3 vials, freeze - dried
Storage conditions	+2 to +8 °C
Test Conditions	The prepared mycoplasma particles of 2 vials are rehydrated in 1 ml of cell culture medium (DMEM with 5 % FCS) and 10 mM Tris.HCl pH 8.5. Both samples are DNA extracted using the Microsart Extraction or Venor®GeM Sample Preparation Kit. The samples are used as template for mycoplasma PCR with the Microsart ATMP Mycoplasma Kit for qPCR.
Acceptance Criteria, Certified Values	All preparations must show detectable amplification within the defined sensitivity limits. Microsart ATMP Mycoplasma resulted in an average Ct value of 35.3

The manufactured product fulfils all accepted criteria. The lot was released for sale. This certificate is a declaration of analysis at the time of manufacture. This Certificate of Analysis does not exempt the customer from incoming goods inspection.

Certified by Natalie Alkan
Signature:

