
CERTIFICATE OF ANALYSIS

Product Mycoplasma Off® - Mycoplasma Decontamination Spray

Cat. No. 15-1000/15-5000

Lot No. **150M1047**

Date of Expiry 31.10.2018

Date of Analysis 26.04.2017

Quantity 1 liter spray bottle / 5 x 1 liter refill bottle

Store at ambient temperature

Concentration of main active component:

- Glutaral 0.11 [%]

Chemical / physical analysis:

Results of Analysis

- specific weight: 0.904 [g/ml]
- pH value: 7.7
- refraction index: 1.3678
- turbidity: 0.3 [TE/F]

The manufactured product fulfils all accepted criteria. The lot was released for sale. This certificate is a declaration of analysis at the time of manufacture. This Certificate of Analysis does not exempt the customer from incoming goods inspection.

Certified by Natalie Alkan
Signature:

